
	<b>DEPARTMENT of PARASITOLOGY</b> <b>BPRC Rijswijk</b>	SOP: 001 Version: 001
	<b>Standard Operating Procedure for</b> <b>AMA1, HBsAg and Ag85A formulations in Alhydrogel®, Squalene-in-</b> <b>water and QS21-Liposomes</b>	Page 1 of 4  Date: 19/10/12
<b>Written by: Sumera Y. Younis, MSc</b> <b>Reviewed by: Ed J. Remarque, PhD</b>		

<b>Index</b>	<b>Page</b>
<b>A.1. Aim</b>	<b>2</b>
<b>A.2. Brief description</b>	<b>2</b>
<b>B.1. Equipment and instrument</b>	<b>2</b>
<b>B.2. Materials</b>	<b>2</b>
<b>C.1. Preparation adjuvants</b>	<b>2</b>
<b>C.2. Preparation antigens</b>	<b>2</b>
<b>D.1. Formulations for control vaccines</b>	<b>3</b>
<b>D.2. Formulations for mouse experiments</b>	<b>3</b>

ORIGINAL

	<b>DEPARTMENT of PARASITOLOGY</b> <b>BPRC Rijswijk</b>	SOP: 001 Version: 001
	<b>Standard Operating Procedure for</b> <b>AMA1, HBsAg and Ag85A formulations in Alhydrogel®, Squalene-in-</b> <b>water and QS21-Liposomes</b>	Page 2 of 4
	<b>Written by: Sumera Y. Younis, MSc</b> <b>Reviewed by: Ed J. Remarque, PhD</b>	Date: 19/10/12

### **A.1. Aim**

The formulation of AMA1, HBsAg and Ag85A in Alhydrogel®, Squalene-in-water and QS21-Liposomes

### **A.2. Brief description**

The formulations described here yield 20 µg for AMA1, 100 µg for HBsAg and 200 µg for Ag85A at an injection volume of 1 mL. Per immunisation 50 µL injections were prepared.

### **B.1. Equipment and instruments**

- ❑ Laminar Flow Cabinet: Clean Air
- ❑ Terumo Insulin syringes 0.5 mL: Multimed # 14.BS=05M2713
- ❑ 15 mL tubes: Greiner Bio-one # 188271

### **B.2. Materials**

- ❑ AMA1 vials (62.5 µg per vial; 3 mL vials): BPRC
- ❑ Ag85A vial (228 µL; 2192 µg/mL per vial): Lionex GmbH # LRP-0004.4
- ❑ HBsAg vial (6.56 mg/mL per vial; 100 µL): VFL, Lausanne
- ❑ NaCl 0.9% (Saline for injection): 9501C12  
Use unopened bottle, discard after use.
- ❑ Alhydrogel® 2% (AH 1173 µg; 0.345 mL): Brentag Biosektor
- ❑ SWE (SWE LE 009; 0.345 mL): VFL, Lausanne
- ❑ Cholesterol-DOPC Liposomes (LIPO-003; 0.138 mL): VFL, Lausanne
- ❑ LS saponin (LS-001; 0.138 mL): VFL, Lausanne

### **C.1. Preparation adjuvants**

#### *QS21-liposomes preparation*

1. Take 135 µL cholesterol-DOPC liposomes
2. Add 135 µL of LS saponin
3. Add 67.5 µL saline and mix gently

*AlOH and SWE are ready-to-use*

### **C.2. Preparation antigens**

#### *AMA1 dilution @ 40 µg/mL*


1. Reconstitute one vial containing 62.5 µg in 1562.5 µL saline.
2. Final concentration: 40 µg/mL

#### *HBsAg dilution @ 200 µg/mL*

1. Take 25.9 µL HBsAg @ 6.56 mg/mL
2. Add 824.1 µL saline and mix gently
3. Final concentration: 200 µg/mL

#### *Ag85A dilution @ 400 µg/mL*

1. Take 1 vial of Ag85A and take out 200 µL and transfer to 1.25 mL vial
2. Add 896 µL saline and mix gently

	<b>DEPARTMENT of PARASITOLOGY</b> <b>BPRC Rijswijk</b>	SOP: 001 Version: 001
	<b>Standard Operating Procedure for</b> <b>AMA1, HBsAg and Ag85A formulations in Alhydrogel®, Squalene-in-</b> <b>water and QS21-Liposomes</b>	Page 3 of 4
	<b>Written by: Sumera Y. Younis, MSc</b> <b>Reviewed by: Ed J. Remarque, PhD</b>	Date: 19/10/12

3. Final concentration: 400 µg/mL

### **D.1. Formulations for control vaccines**

#### *AIOH control*

1. Take 40 µL AIOH.
2. Add 40 µL saline and mix gently.
3. Incubate for 1 hour at RT.
4. Take up 50 µL of the mixture into a labeled insulin syringe for one mouse.

#### *SWE control*

1. Take 40 µL SWE.
2. Add 40 µL saline and mix gently.
3. Take up 50 µL of the mixture into a labeled insulin syringe for one mouse.

#### *QS21-Liposomes control*

1. Take 40 µL QS21-liposomes.
2. Add 40 µL saline and mix gently.
3. Take up 50 µL of the mixture into a labeled insulin syringe for one mouse.

### **D.2. Formulations for mouse experiments**

#### **AMA1**

##### *Adjuvant-Formulations AMA1 - AIOH*

1. Take 270 µL AMA1 at 40 µg/mL.
2. Add 270 µL AIOH and mix gently.
3. Incubate for 1 hour at RT.
4. Fill labeled insulin syringes with 50 µL of the mixture for each mouse.

##### *Adjuvant-Formulations AMA1 - SWE*

1. Take 270 µL AMA1 at 40 µg/mL.
2. Add 270 µL SWE and mix gently.
3. Fill labeled insulin syringes with 50 µL of the mixture for each mouse.


##### *Adjuvant-Formulations AMA1 - QS21-Liposomes*

1. Take 270 µL AMA1 at 40 µg/mL.
2. Add 270 µL QS21-liposomes and mix gently.
3. Fill labeled insulin syringes with 50 µL of the mixture for each mouse.

#### **HBsAg**

##### *Adjuvant-Formulations HBsAg - AIOH*

1. Take 270 µL HBsAg at 200 µg/mL.
2. Add 270 µL AIOH and mix gently.
3. Incubate for 1 hour at RT.
4. Fill labeled insulin syringes with 50 µL of the mixture for each mouse.

	<b>DEPARTMENT of PARASITOLOGY</b> <b>BPRC Rijswijk</b>	SOP: 001 Version: 001
	<b>Standard Operating Procedure for</b> <b>AMA1, HBsAg and Ag85A formulations in Alhydrogel®, Squalene-in-</b> <b>water and QS21-Liposomes</b>	Page 4 of 4
	<b>Written by: Sumera Y. Younis, MSc</b> <b>Reviewed by: Ed J. Remarque, PhD</b>	Date: 19/10/12

*Adjuvant-Formulations HBsAg - SWE*

1. Take 270  $\mu\text{L}$  HBsAg at 200  $\mu\text{g}/\text{mL}$ .
2. Add 270  $\mu\text{L}$  SWE and mix gently.
3. Fill labeled insulin syringes with 50  $\mu\text{L}$  of the mixture for each mouse.

*Adjuvant-Formulations HBsAg - QS21-Liposomes*

1. Take 270  $\mu\text{L}$  HBsAg at 200  $\mu\text{g}/\text{mL}$ .
2. Add 270  $\mu\text{L}$  QS21-liposomes and mix gently.
3. Fill labeled insulin syringes with 50  $\mu\text{L}$  of the mixture for each mouse.

Ag85A

*Adjuvant-Formulations Ag85A - AIOH*

1. Take 270  $\mu\text{L}$  Ag85A at 400  $\mu\text{g}/\text{mL}$ .
2. Add 270  $\mu\text{L}$  AIOH and mix gently.
3. Incubate for 1 hour at RT.
4. Fill labeled insulin syringes with 50  $\mu\text{L}$  of the mixture for each mouse.

*Adjuvant-Formulations Ag85A - SWE*

1. Take 270  $\mu\text{L}$  Ag85A at 400  $\mu\text{g}/\text{mL}$ .
2. Add 270  $\mu\text{L}$  SWE and mix gently.
3. Fill labeled insulin syringes with 50  $\mu\text{L}$  of the mixture for each mouse.

*Adjuvant-Formulations Ag85A - QS21-Liposomes*

1. Take 270  $\mu\text{L}$  Ag85A at 400  $\mu\text{g}/\text{mL}$ .
2. Add 270  $\mu\text{L}$  QS21-liposomes and mix gently.
3. Fill labeled insulin syringes with 50  $\mu\text{L}$  of the mixture for each mouse.

ORIGINAL